




TAGMe DNA methylation detection

Female Genital Tract Cancer Theranostics
Urinary Tract Cancer Theranostics

EPIPROBE

>>> Why receive methylation detection?





Dr. Wenqiang Yu

- Doctoral Supervisor of Fudan University, Chief Scientist of national "973" project,
- Changjiang Scholar Distinguished Professor, PI of Epigenetics Center of Fudan Biomedical Research Institute.
- 2001-2007 Postdoctoral Fellow, Uppsala University, Sweden; **Johns Hopkins University, USA (tutored by Dr. Andy);**
- In November 2007, Faculty and Associate Research Scientist of Columbia University.



Research findings


The whole genome DNA methylation sequencing method with independent intellectual property right, **Guided Positioning Sequencing (GPS) technology**, was established and **Tumor Aligned General Methylated Epiprobe (TAGMe)** were discovered, which has been double-blind verified in **50000+** clinical samples.

>>> Discovery of TAGMe

Roadmap
2003~2012

ENCODE
2003~2012

TCGA
2006~2018




Methylation
Calculation



Guide
Positioning



Research: From simple to complex




Application: From complex to simple



GPS method: covering 96% of the genome C (1.12G/1.17G);
WGBS method: covering 60 to 80% CpG of the genome.




Tumor Aligned General Methylated Epiprobe (TAGMe)



Genome Research, 2019.01 (**IF=11.1**)
Cancer Research, 2019.10 (**IF=12.7**)
Clinical and Translational Medicine, 2021.06 (**IF=11.5**)
Frontiers in Molecular Biosciences (**IF=5.2**)
Signal Transduction and Targeted Therapy (**IF=38.104**)

>>> Technological breakthrough: TAGMe-DNA methylation detection method

Establish a whole-process standardized system from sample collection to report generation



Me-qPCR platform

Without bisulfite treatment
Automate DNA methylation
detection in one step

Applicable detection items


- Cervical cancer,
- Urinary tract cancer,
- Endometrial cancer

Technical advantages without bisulfite treatment


- More stable
- More sensitive
- More convenient
- Automated

Original DNA methylation detection technology, based on qPCR platform that doesn't require bisulfite treatment
-- Me-qPCR, can detect as low as 0.2% of tumor components.

Me-qPCR vs Pyro-Seq




Limit of Detection: 0.2%



>>> Double-blind validation of TAGMe

No.	Cancer type	Sample type	Sample Number	Partial research results	Performance
1	Cervical cancer	Cervical cells/vaginal secretions, etc.	36348	Clin and Trans Med (2021.06, the latest IF is 11.5) STTT (2022.07, IF 38.1)	94.3% specificity, 96% sensitivity (exfoliate cells)
2	Urinary tract carcinoma	Urine/tissue, etc.	3499	Article is scheduled for publication in 2022	92.7% specificity, 82.1% sensitivity (urine)
3	Lung cancer	Alveolar lavage fluid/pleural fluid/tissue/blood, etc.	3385	Cancer Res (2019.10, the latest IF is 12.7)	96.5% specificity, 87% sensitivity (lavage fluid/pleural fluid)
4	Endometrial cancer	Cervical/uterine cavity cells/tissues, etc.	884	Front Mol Biosci (2021.11, the latest IF is 5.2)	87.3% specificity, 90.9% sensitivity (exfoliate cells)
5	Biliary tract tumors	Bile/tissue, etc.	930	Article is scheduled for publication in 2022	100% specificity, 96.9% sensitivity (tissue)
6	Immunotherapy	blood	746	Article is scheduled for publication in 2022	85.7% specificity, 66.7% sensitivity (Blood)
7	breast cancer	Tissue/blood, etc.	150	R&D is in progress	92.5% specificity, 100% sensitivity (tissue)
8	Liver cancer	Tissue/blood, etc.	979	Article is scheduled for publication in 2022	90.1% specificity, 82.2% sensitivity (tissue)
9	gastric cancer	Tissue/blood, etc.	196	R&D is in progress	100% specificity, 90% sensitivity (tissue)
10	Colorectal cancer	Feces/tissue/blood, etc.	189	R&D is in progress	90% specificity, 100% sensitivity (tissue)
11	Thyroid cancer	Tissue	215	R&D is in progress	-
12	Other	Cerebrospinal fluid/bone marrow smear/tissue/blood, etc.	971	R&D is in progress	-

Summary: By the end of March 2022, total of 50552 clinical samples have been double-blind validated, and the overall consistency rate of tissue samples is >90%.






TAGMe
DNA methylation
detection for female
genital tract cancer



Cervical Cancer(TAGMe-CeCan) & Endometrial
Cancer(TAGMe-EnCan)

Eliminate the cancer in the precancerous stage

- According to the analysis of the 3728 newly enrolled and unblinded samples in the double-blind verification, performance of DNA methylation detection is as follows:
pathologic positive includes high-grade cervical lesions and cervical cancer patients, and pathology-negative refers to diseases that have not reached high-grade lesions or cervical cancer.




Clinical and Translational Medicine.


2021.06 (IF=11.5)

- DNA methylation in cervical high-grade lesions (HSIL) screening has a specificity of **95.2%**, and sensitivity of **76.5%**;
- Overall (\geq HSIL) screening specificity is **95.2%**, and sensitivity is **84.4%**.


Results of double-blind validation

Further analysis of the detection performance of DNA methylation detection in different stages found that the detection rate in the cancer group was 100%, and the detection rate in the high-grade lesion group was 76.5%.





Early endometrial cancer is confined to the uterus



Endometrial biopsy may result in discomfort, bleeding, infection, and uterine perforation, with a high rate of missed tests

- Pain point: Lacking sensitive and accurate non-invasive screening method. Symptoms such as early irregular vaginal bleeding and vaginal drainage are easily overlooked, missing the opportunity for early diagnosis.

■ Transvaginal ultrasonography:

Convenient and noninvasive, vaginal ultrasound is easy to miss diagnosis when the endometrium is <5 mm thick and hard to assess premenopausal endometrial lesions.

■ Hysteroscopy:

Expensive, most patients require anesthesia, and has the risk of side effect(e.g. infection, water intoxication, air embolism, etc.), thus cannot be applied as a routine screening method.

■ Microscopic diagnostic curettage:

Invasive surgery, significant pain, has the risk of side effect(e.g. bleeding, infection, uterine perforation, uterine adhesions etc.), along with the possibility of missed scratch.

■ Endometrial biopsy:

Gold standard, invasive surgery, and it is easy to get insufficient and inaccurate samples, especially for postmenopausal patients



>>> Product performance: endometrial cancer

Double-blind samples:

- ① Pap sample (109) : NE 47 , EH 20 , AH 9 , EC 33 .
- ② Tao sample (103) : NE 44 , EH 20 , AH 8 , EC 31 .

Diagnostic model

Diagnostic model	Endometrial cancer detection performance				
	Sensitivity (%)	Specificity (%)	PPV (%)	NPV (%)	Accuracy (%)
Pap Brush (Cutoff:<4.03)	80.65	82.81	86.49	62.69	84.21
Tao Brush (Cutoff:<1.25)	61.29	95.31	54.05	87.50	82.11
Either positive as positive (Cutoff: Pap<4.03, Tao<1.32)	90.32	73.44	62.22	94.00	78.95
Both positive as positive (Cutoff: Pap<2.5, Tao<4.55)	61.29	100.00	100.00	84.21	87.37




NE: Normal endometrium

EH: Endometrial hyperplasia

AH: Atypical hyperplasia

EC: Endometrial cancer




Combined with the existing clinical uterine cavity and cervical exfoliation cell sampling devices - Pap brush and Tao brush, the non-invasive screening of endometrial cancer is realized.


Double-blind test results show:

Pap brush: AUC =0.86, specificity=82.81%, sensitivity=80.65%;

Tao brush: AUC=0.83, specificity=95.31%, sensitivity=61.29%;

Project	Application scenarios	Sample types	Sample volume	Transportation requirement
TAGMe DNA methylation detection for female genital tract cancer -	<ul style="list-style-type: none"> • Early cancer screening • Auxiliary diagnosis • Risk monitoring • Recrudescence monitoring 	Cervical scraping (TCT sampling method)	2~5 mL	Room temperature






Full-process solution for urinary tract cancer

TAGMe-UrCan





Urine detection for bladder cancer



Double-blind validation of urine samples for bladder cancer,
AUC=0.93, specificity = 92.9%, sensitivity = 85.7%;




- Bladder cancer: 90% of which is transitional epithelial cell carcinoma; 80% of which is non-invasive carcinoma; 93% of which is low-grade carcinoma.



- The detection rate of high-grade bladder cancer is **92.8%** and that of low-grade bladder cancer is **82.7%**;
- The detection rate of invasive carcinoma (MIBC) is **95.9%**, The detection rate for non-invasive carcinoma (NMIBC) is **87.8%**.

Urine detection for renal pelvis cancer and ureteral cancer



- Renal pelvis cancer: sensitivity=91.7%, specificity=92.1%;
- Ureteral cancer: sensitivity=78.7%, specificity=92.1%;




>>> Application scenario: Auxiliary diagnosis

Application scenario 1: Auxiliary diagnosis is carried out by TAGMe-UrCan for patients with suspected bladder cancer in outpatient clinics, inability to do cystoscopy due to urethral stenosis, suspected upper urothelial carcinoma, etc.;

Positive: high risk of cancer, recommend to receive cystoscopy to confirm, or follow the doctor's advice;

Negative: low risk of cancer, regular test of DNA methylation tests, or follow the doctor's advice;



Application scenario 2: For patients with clinical symptoms such as space-occupying lesions, hematuria, etc., but the cystoscopy comes back negative;


Positive: high risk of cancer, recommend to re-examine for cystoscopy, close monitoring, and follow the doctor's advice;

Negative: low risk of cancer, regular check-up, cancer risk monitoring;




Case study of surgical efficacy determination

All patients (12/12) had lower postoperative DNA methylation level than before, but postoperative DNA methylation level of yellow SS was still a strong positive, suggesting a high suspicion of postoperative tumor residue. It was strongly recommended to review and monitor.



Application Scenario 3: Evaluation of Surgical Efficacy:

1. Methylation value is still above the threshold (positive) after the surgery: recommend to carry out perfusion chemotherapy and closely monitor the risk of recrudescence;
2. Methylation value is below the threshold (negative) after the surgery: follow the doctor's comprehensive judgements for follow-up treatment plan;



Application scenario 4: Evaluate the efficacy by dynamic changes in methylation values after each perfusion chemotherapy

1. The postoperative methylation value is still above the threshold (positive): recommend to continue (or change the regimen) for perfusion chemotherapy and closely monitor the risk of recrudescence;
2. The postoperative methylation value is below the threshold (negative): follow the doctor comprehensive judgements on whether to stop perfusion chemotherapy;

Evaluate the efficacy based on changes in methylation values



>>> Recrudesce of urinary urothelial carcinoma

Case studies of recrudesce monitoring

Clinical procedures

Patient: Lu **, male, 82 years old
First diagnosis: malignancy of the bladder

21.01.04
Cystoscopy: high-grade invasive urothelial carcinoma with squamous metaplasia

21.01.19
Surgery: Transurethral resection for bladder cancer (TURBT)

21.05.10
Cystoscopy: 1.5*0.6 cauliflower-like neoplasia is seen at the lateral margin of the right ureteral orifice, indicating recrudesce

Outpatient service

Diagnosed by cystoscopy

Surgical treatment

Leave hospital

TAGMe detection-1

TAGMe detection-2

TAGMe detection-3

Surgical treatment

TAGMe detection-1

Continue follow-up

21.01.18
Before surgery:
TAGMe Value = 0.1, **Positive**

21.02.18
TAGMe Value = 3.5, negative, indicating no recrudesce


21.03.16
TAGMe Value = 1.8, **positive, indicating recrudesce**

21.04.19
TAGMe Value = -0.2, **positive, indicating recrudesce**

21.05.26
TAGMe Value = 4.7, negative, indicating no recrudesce

Recrudesce monitoring

TAGMe Methylation detection



Application scenario 5: Within 3 months after surgery, receive TAGMe-UrCan once a month can dynamically detect the risk of early postoperative recrudescence

1. Positive: high risk of recrudescence, recommend to have cystoscopy review early, and closely monitor;
2. Negative: the risk of recrudescence is low, follow the doctor's comprehensive judgment and advice.

Application scenario 6: After 3 months of surgery, before each cystoscopy, have TAGMe-UrCan to detect the risk of recrudescence after surgery for long-term dynamic detection

1. Positive: high risk of recrudescence, recommend to have cystoscopy review early, and closely monitor;
2. Negative: the risk of recrudescence is low, follow the doctor's comprehensive judgment and advice.



>>> Changhai Hospital - Full Process Case Study

Clinical procedures

Patient: Liu **, female, 82 years old
Initial diagnosis: high-grade varus papillary urinary tract carcinoma of the bladder

20.03.31 Cystoscopy: Cauliflower-like neoplasia is seen at the ureteral opening on the left side of the bladder, and cauliflower new organism is seen at the top of the bladder

20.05.22 Surgery: Transurethral resection for bladder cancer (TURBT)

20.06~20.12: Seven months after the operation, due to the inconvenience of her location, cystoscopy pain and other reasons, the patient's compliance was poor, and she did not return to the hospital for re-examination



21.01.27 Cystoscopy: Cauliflower-like neoplasia on the anterior wall of the bladder; 21.02.01 Pathology: High-grade papillary urothelial carcinoma

Outpatient service

Diagnosed by cystoscopy

Surgical treatment

Leave hospital

TAGMe recrudescence detection-1

TAGMe recrudescence detection-2

Cystoscopic review

20.03.24
TAGMe Value = 92, **positive**, indicating a high risk of urinary urothelial cancer

20.05.22
Before surgery:
TAGMe Value = 94, **positive**

20.05.22
After surgery:
TAGMe Value = 42, negative, indicating successful surgery

20.12.26
TAGMe Value = 85, **positive**, indicating recrudescence

21.01.26
TAGMe Value = 95, **positive**, indicating recrudescence and progression; doctors strongly recommend returning to the hospital for re-examination!

Auxiliary diagnosis

Efficacy assessment

recrudescence monitoring

TAGMe-UrCan

- Sampling tubes are mailed directly to patients. Urine can be taken at home, which is convenient!
- TAGMe can be applied to whole process from early screening, auxiliary diagnosis, efficacy evaluation to recrudescence monitoring of urinary urothelial cancer.





Sample requirements

Noninvasive: Only 30 ml of urine is required

Convenient: Samples can be taken at home

Simple: Storage and shipping at room temperature for 15 days

Fast: Electronic reports takes only 3 to 5 work days

Scenarios

High-risk populations, early screening

Clinical-suspicious people, assist cystoscopy for diagnosis


Perfusion chemotherapy, efficacy assesment

Initial prognosis, risk assessment of recrudescence

Prognostic recrudescence, long-term recrudescence monitoring

Epiprobe TAGMe DNA methylation molecular detection provides a full-process solution for urinary tract carcinoma





Before diagnosis:-Adjunctive diagnosis and recurrence monitoring

During treatment. : Evaluation of surgical efficacy and perfusion therapy efficacy

After surgery and long-term recurrence monitoring





Genome Research, 2019.01 (**IF=11.1**)
Cancer Research, 2019.10 (**IF=12.7**)
Clinical and Translational Medicine, 2021.06 (**IF=11.5**)
Frontiers in Molecular Biosciences (**IF=5.2**)
Signal Transduction and Targeted Therapy (**IF=38.104**)



About Epiprobe

As a high-tech enterprise founded in 2018 by top epigenetic experts, Epiprobe focuses on the molecular diagnosis of cancer DNA methylation and precision theranostics industry. With a profound technology basis, EPIPROBE aims to lead the era of new products to nip cancer in the bud!

Based on Epiprobe core team's long-term research, development and transformation in the field of DNA methylation with the cutting-edge innovations, combined with the unique DNA methylation targets of tumors, Epiprobe uses a unique multivariate algorithm combining big data and artificial intelligence technology to independently develop an exclusive patent-protected liquid biopsy technology. By analyzing the methylation level of specific sites of free DNA fragments in the sample, the shortcomings of traditional examination methods and the limitations of surgery and puncture sampling are avoided, which not only achieves accurate detection of early tumors, but also enables real-time monitoring of tumor occurrence and development dynamics.

Epiprobe's tumor molecular detection technology can be used for early tumor screening, auxiliary diagnosis, preoperative and postoperative evaluation, recrudescence monitoring, which runs through the whole process of tumor diagnosis and treatment, providing better solutions for doctors and patients.

Add: 6th Floor, Building 4 , No. 333 Guiping Road, Xuhui District, Shanghai,China
Email: service@epiprobe.com

